



The **Philips Heartstart MRx** is a lightweight, portable monitor/defibrillator. with the capability to offer four different modes of operation; monitor, manual Defibrillate, AED, and pacing. The defibrillator was designed to provide both monitoring and resuscitation by using mutual-parameter functions, including algorithms and options for transferring data to, a full range of defibrillation therapies. The Heartstart MRx is fast and reliable with 12-lead ECG and advanced STEMI support. The MRx has a large 8-inch color screen that can display up to 4 waveforms and numerics. The Philips Heartstar MRx additionally monitors other parameters: arrhythmia detection, SpO2, NIBP, IBP, EtCO2, and temperature.

Features

- Easy to use Interface with a large color display
- Can Display 4 waveforms and numerics, or all 12 leads and numerics
- Screen views can be customized to fit the user's preferences
- Normal or high-contrast view of easy viewing in different lighting conditions
- Seamless data transmission for 12 lead ECGs with optimal wireless link
- 10 hours of monitoring with 2 fully charged batteries
- Ready-For-Use indicator
- optimal AC or DC battery options



Specifications

Dimensions

Without External Paddles

Width: 12.4 in (313 mm)
Depth: 8.3 in (210 mm)
Height: 11.7 in (295 mm)

With External Paddles

Width: 13.4 in (340 mm)
Depth: 8.3 in (210 mm)
Height: 13.6 in (345 mm)

Weight

13.2 lbs. (6kg): Base unit with 1 battery, pads, and pads cable.

+4.1 lbs. (1.86 kg): With added carrying case.

+2.5 lbs. (1.1 kg): With added paddle tray and external standard paddles.

Environmental

Water Resistance: Meets IEC 60601-2-4

Solids Resistance: IP2X

Temperature

Operating: 32° - 113° F (0° - 45° C)

Storage: -4° - 158° F (-20° - 70° C)

Humidity: Operating: 0% to 95% relative

Safety: Meets EN 60601-1, UL 2601-1, CSA C22.2 No. 601-1-M90 CSA, EN 60601-2-4

Display

Dimensions: 8.4 in diagonal (128 mm x 171 mm)

Type: TFT color LCD

Resolution: 640 x 480 pixels (VGA)

Wave Viewing Time: 5 seconds (ECG)

Strip Chart Printer

Printer

Standard: 50 mm (paper width) thermal array printer.

Optional: 75 mm (paper width) thermal array printer.

Continuous ECG Strip: Prints primary ECG lead with event annotations and measurements in real-time or with 10-second delay.

Auto Printing: Printer can be configured to print marked events, charge, shock, and alarms.

Reports: Event Summary, 12-Lead, Vital Signs Trending, Operational Check, Configuration, Status Log, and Device Information.

Paper Size

1.97 in (50 mm) W by 100 ft. (30 m) L.

2.95 in (75 mm) W by 100 ft. (30 m) L.



Specifications Continued

Defibrillator

Model: HeartStart MRx (M3535A)

Waveform: Biphasic Truncated Exponential. Waveform parameters adjusted as a function of patient impedance.

Output Energy: Manual (selected): 1-10, 15, 20, 30, 50, 70, 100, 120, 150, 170, 200 Joules maximum energy, limited to 50 Joules for internal defibrillation. AED Mode (single energy output): 150 Joules into a 50 ohm load.

Charge Time: Less than 5 seconds to 200 Joules with a new, fully charged lithium ion battery at 25° C.

Shock Delivery: Via multifunction defib electrode pads or paddles.

Quick Shock: Less than 10 seconds from cessation of CPR to shock delivery.

Patient Impedance Range

Minimum: 15 ohm (internal defibrillation); 25 ohm (external defibrillation).
Maximum: 150 ohm.

AED Mode: Shock advisory sensitivity and specificity meet AAMI DF-39 guidelines.

Battery

Type: 6.0 Ah, 14.8 V, rechargeable lithium ion.

Dimensions

Height: 6.5 in (165 mm)

Width: 3.8 in (95 mm)

Depth: 1.6 in (42 mm)

Weight: 1.6 lbs. (0.73 kg)

Charge Time: Approximately 3 hours to 100%, 2 hours to 80%.

Capacity

At least 5 hours of monitoring with ECG, SpO₂, CO₂, temperature and two invasive pressures monitored continuously, NBP measured every 15 minutes, and 20 200J discharges (with a new, fully charged battery, operating at room temperature, 25° C).

At least 3.5 hours of monitoring with ECG, SpO₂, CO₂, temperature and two invasive pressures monitored continuously, NBP measured every 15 minutes, and pacing at 180ppm at 160mA.

Battery Indicators

Battery gauge on battery, capacity indicator on display; flashing RFU indicator, chirp, and 'Low Battery' message appears on display for low battery condition, when 10 minutes of monitoring time and 6 maximum energy discharges remain (with a new battery at room temperature, 25° C).



Specifications Continued

Data Storage

Internal: 12 hours of continuous ECG waveforms and events, maximum capacity of 55 event summaries.

Data Card: 60 event summary reports or 240 megabytes of patient data.

ECG and Arrhythmia Monitoring

Input

Up to 4 ECG waves displayed and up to 2 ECG waves print simultaneously.

Lead I, II, or III obtained through 3-lead ECG cable and separate monitoring electrodes. With 5-lead cable, obtain leads aVR, aVL, aVF, or V. Pads ECG obtained through 2 multifunction defibrillation electrode pads.

Lead Fault: 'Lead Off' message and dashed line displayed, if an electrode or lead wire becomes disconnected.

Pads Fault: Dashed line displayed if a pad becomes disconnected.

Heart Rate Display: Digital readout on display 15 to 300 bpm, accuracy $\pm 10\%$.

Heart Rate . Arrhythmia Alarms: HR, Asystole, VFIB/VTACH, VTACH, extreme tachycardia, extreme bradycardia, PVC rate, Pacer not capture, Pacer not pacing.

ECG Size: 2.5, 5, 10, 20, 40 mm/mV, autogain.

Available Options

Non-Invasive Pacing: SpO2 pulse oximetry.

Non-Invasive Blood pressure: CO2 monitoring.

Invasive Blood Pressure (2 lines): Continuous temperature monitoring.

12-Lead Acquisition: 12-lead transmission.

Q-CPR Measurement and Feedback: Audio recording.

ACI-TIPI & TPI Predictive Instruments: Periodic clinical data transmission.

Batch/LAN Data Transfer

