



The **Philips Respironics Trilogy 202** is a Ventilator that offers both invasive and non-invasive ventilation. The Trilogy ventilator can control volume and pressure in the ventilation, making it easier for the patient to breath. The Trilogy is a ventilator that has a microprocessor control built into the system. Unlike the 100 and 200 Version of the Trilogy, it has the air and Oxygen blended in the machine and runs at 50 psi, where the other version of the Philips Respironics Trilogy run off of low flow O2. The machine can be used in both home care, hospital care, and hospital transport settings. It comes with an internal battery that can last up to 3 hours and an additional backup battery that lasts up t 6 hours. The Philips Respironics Trilogy 202 can be used on a wide range of patients from Pediatrics to Adults.

The Respironics Trilogy provides continuous ventilation support for those individuals that are suffering from acute or chronic respiratory failure, and chronic repository distress. The system comes with an Auto-Track algorithm that can enhance the patients breathing by adjusting and changing the breathing patterns in the ventilator.

The Philips Respironics Trilogy ventilation system comes with an SD card slot on the side of the machine. it can record breath by breath patient data and the last 72 hours of waveforms This makes it easy and convenient to transfer the Data from the ventilator to a central data station in the hospital setting.

## Features

- Compact design with long-life internal battery for intra-hospital transport.
- Supports active and passive circuits for invasive and noninvasive ventilation.
- Innovative leak compensated volume control ventilation.



# Specifications

## Controls

CPAP : 4 - 20 cmH2O (passive leak port circuits) cm H2O  
 IPAP: 50 cm H2O  
 EPAP/PEEP (a): 0 - 25 cmH2O (active valve circuits) cm H2O  
 EPAP/PEEP (b): 4 - 25 cmH2O (passive leak port circuits) cm H2O  
 Pressure support: 0 – 30 cmH2O cm H2O  
 Tidal volume: 50 - 2000 ml  
 Breath rate(a): 0 - 60 BPM (AC mode) beats per minute  
 Breath rate(b): 1 - 60 BPM (all other modes) beats per minute  
 Inspiratory time: 0.3 - 5.0 s Synchrony features  
 Rise time: 1 - 6 (relative scale)  
 Ramp start pressure (a): 0 - 25 cmH2O (active circuits) cm H2O  
 Ramp start pressure (b): 4 - 25 cmH2O (passive circuits) cm H2O  
 Ramp start pressure (c): 4 - 19 cmH2O (CPAP mode) cm H2O  
 Ramp time Off: 5 - 45 min  
 C-Flex Off: 1 - 3 (relative scale)  
 Flow trigger sensitivity: 1 - 9 l/min  
 Flow cycle sensitivity: 10 - 90 %

## Environmental

Operating temperature: 5 – 40 °C  
 Relative humidity: 15 - 95 %  
 Storage temperature: -20 – 60 °C  
 Atmospheric pressure: 60 - 110 kPa (450 - 825 mmHg)

## Electrical

Input voltage: 100 - 240 VAC and 50/60 Hz and 2.1 A  
 Detachable battery voltage: 14.4 VDC  
 Internal battery life: 3 hours under normal conditions hr  
 Detachable battery life: 3 hours under normal conditions hr  
 external battery connection: 12 VDC

## Volume Modes

MPV Mouthpiece Ventilation Mode: Yes  
 Pressure control (PC): Yes  
 Synchronized intermittent mandatory ventilation (SIMV): Yes  
 SIMV with pressure support (SIMV w/PS): Yes  
 Control ventilation (CV): Yes

## Alarms

Low tidal volume Off: 50 - 2000 ml  
 Circuit disconnect Off: 10 - 60 sec  
 Apnea Off: 10 - 60 s and 4 - 60 BPM  
 High tidal volume Off: 50 - 2000 ml  
 High minute ventilation Off: 1 - 99 l/min Off, 1 - 99 l/min

# Specifications

## Alarms (continued)

Low minute ventilation Off: 1 - 99 l/min  
 High respiratory rate Off: 4 - 80 beats per minute  
 Low respiratory rate Off: 4 - 80 beats per minute

## Pressure Modes

Pressure control (PC): Yes  
 Pressure control-SIMV (PC-SIMV): Yes  
 Spontaneous ventilation (S): Yes  
 Spontaneous ventilation with timed back-up (S/T): Yes  
 Timed ventilation (T): Yes

## Options

Detachable back-up battery: Offering up to 3 hours additional operating time  
 Hospital roll stand: Provides convenient accessory basket and humidifier mount  
 DirectView: 1 GB SD card data storage integrated into the ventilator software

## Monitored Parameters

Tidal volume: 0 - 2000 ml  
 Minute ventilation: 0 - 99 l/min  
 Estimated leak rate: 0 - 200 l/min  
 Respiratory rate: 0 - 80 beats per minute  
 Peak inspiratory flow: 0 - 200 l/min  
 Peak inspiratory pressure: 0 - 99 cm H<sub>2</sub>O  
 Mean airway pressure: 0 - 99 cm H<sub>2</sub>O  
 % patient triggered breaths: 0 - 100 %  
 I:E ratio: 9.9:1 - 1:9.9

## Synchrony Features

Auto-Trak sensitivity: Auto-Trak sensitivity Auto-adaptive triggering, cycling, and leak compensation (available in all modes, passive circuit only)  
 Adjustable flow triggering: 1 – 9 l/min (available in all modes and circuit types) l/min

## Circuit Types

Active exhalation valve with proximal airway pressure (PAP): Yes  
 Active exhalation valve with flow sensor: Yes  
 Passive exhalation port: Yes

