

Medfusion 3010a Syringe Pump



Typical Manufacturer's Picture

- 16 infusion modes.
- Graphic pressure trend display.
- Begins infusing in just 3 keystrokes.
- Easy to see screen.
- MRI Compatible.
- User Friendly.
- Stackable horizontal design.

Specifications

Dimensions

Size: 10.5" x 5.1" x 5.5". Weight: 3.8 lbs.

System Accuracy

+/- 2%, excluding syringe variations

Flow Rate

0.01 ml/hr to 1130 ml/hr, dependent upon syringe size

Delivery Modes

ml/hr; Volume / Time; Mass; Body Weight; Custom Dilution; Intermittent

Delivery Options

Volume Limit; Bolus Dose; Loading Dose; KVO Rate; Standby; Override Occlusion Limit

Delivery Libraries

Four libraries with 16 profiles each

Power

ACL 100 – 240 VAC, 50/60 Hz, 30Va
External DC source (optional)

Recharge Time

Less 16 hours with the pump turned off.

Battery Capacity

Rechargeable. 10 hours at 5.0 ml/hr with 60 ml syringe (from fully charged battery)

Soma Technology, Inc. acknowledges all registered trademarks of manufacturers' listed.
The technical data given in this publication are for general information and are subject to change without notice.



166 Highland Park Dr. Bloomfield, CT 06002 U.S.A.

Tel: (860) 218-2575 Fax: (860) 218-2565 Email: soma@somatechnology.com