



Axia Blanket Warmers undergo stringent testing procedures during the manufacturing process to assure long life, the highest quality and accurate temperature settings.

Features

- Keyed temperature lock out
- Push button digital operation
- Large dual digital displays
- Displays actual temperature and set point at all times
- Separate digitally controlled heating chambers
- Rapid warm time (2-6 hours depending on load)
- Uniform heating
- Fully insulated
- Audible and Visual Over-Temp Alarm
- U.L./C.U.L. approved
- FDA Registered
- OSHPD Pre-Appro



Specifications

Overall Size

30"W x 26.5"D x 36"T
2 Adjustable Shelves

Chamber I. D.

26"W x 23"D x 25"T
8.63 cubic feet of storage
Temperature range: 90-160 F.
Accuracy +/- 2 F. degree of set point

Electrical

120V, 7.5 amp., 50/60HZ
.35 kw/hr, 1195 btu/hr
UL/CUL listed
FDA registered

