



The **Axia SWC 72** full size blanket warmer is equipped with push button digital operation, large digital display, and separate digitally controlled heating chambers. This full size blanket warmer is fully insulated and features uniform heating between both chambers. This warming cabinet operates rapidly with a warming time of 2 - 6 hours depending on load.

## Features

- Keypad temperature lock out
- Push button digital operation
- Large dual digital displays
- Displays actual temperature and set point at all times
- Separate digitally controlled heating chambers
- Rapid warm time (2-6 hours depending on load)
- Uniform heating
- Fully insulated
- Audible and Visual Over-Temp Alarm
- U.L./C.U.L. approved
- FDA Registered
- OSHPD Pre-Appro



# Specifications

|                     |  |
|---------------------|--|
| <b>Overall Size</b> | 30"W x 26.5"D x 74.5"T<br>2 Adjustable Shelves<br>2 fixed shelves  |
| <b>Chamber I.D.</b> | 26"W x 23"D x 61"T<br>21 cubic feet of storage<br>Temperature range: 90-160 F.<br>Accuracy: +/- 3 degree F. of set point |
| <b>Power</b>        | 120V, 7.5 amp., 50/60HZ<br>.53 kw/hr, 1810 btu/hr<br>UL/CUL listed<br>FDA registered                                     |

