



The **Cryomedics Circon Cabot MM-4000** is a versatile, high-performance colposcope designed to meet more demanding needs from gynecologic examinations to teaching and research. It features wide field of view binocular eyepieces, an adjustable fine focus, articulating adjustable arm, fiber optic light source, and fiber optic cable.

## Features

- Optical beam splitter for co-observation tube and photographic video.
- Optional 35mm accessories can simultaneously be used for variable magnification and features.
- Fiber-optic light source and cable.
- Wide field of view binocular eyepieces.

