



Model 7924 (glass or solid door) double door warming cabinets (blanket warmers) are 74-1/2" tall and are the perfect answer for warming blankets (blanket warmers) and solutions in two separately heated compartments. Set temperature control via key and pass-code lockout. Self-closing, magnetic gasket, adjustable tension doors. No-set temperature restriction between compartments. Each warming cabinet unit comes standard with integrated casters for ease of mobility.

Features

- Double Compartments
- Comes with Standard Double Glass Door
- Independently Controlled Compartments
- 3 Shelves, open, adjustable design adds convenience and flexibility
- PID Digital controllers feature Tactile Membrane (Stainless Steel Dome)
- Full Numeric Keypads
- Over-Temp audible and Visual Warnings
- Automatic Cut-Off
- Built-In Temperature Data Recorder with Real-Time Clock
- 160 x 180 Pixel Graphical Display with Dual Color (red / green)
- At A Glance: Trend View showing a digital representation of Recent Temperature Conditioning
- Self-Closing, Magnetic Gasket, Adjustable Tension Doors
- Temperature Control Range allows for Setting from 90 to 160 degrees
- Solutions Only warming cabinets can be easily programmed by User to Heat Only Between 90F to 110F
- Self-Tune Temperature feature maintains Set-Temperature to with +/- 1F
- Set temperature can be protected 2 ways: Key lock or User Defined Pass Code
- 2 Standard Door Locks: Key Lock or User Defined Pass Code
- No set Temperature Restriction between compartments
- Calibration Due Date Reminder



Specifications

Dimensions

Height: 74.5in
Width: 30 in
Depth: 26.625 in
Weight: 495 lbs

Shelves

3

Power

Amps:15
Total Wattage: 1500

