



## Features

- Flexible Ureteropyeloscope.
- 160° active deflection, tight bend radius and superior secondary deflection, 9.8 Fr. (3.3 mm)
- AutoSeal™ protection.
- Anti-stretch cable mechanism.
- Access the entire internal collecting system with Anatomy Directed Deflection.
- Accommodates 3 Fr operating instruments, EHL probes, or laser lithotriper fiber, and irrigates at the same time.
- Eliminates vent valves and venting caps for disinfection/sterilization; protects your AUR™ 9 from costly repairs.
- Releases stresses on the deflection cable, maintaining optimum deflection angle over time.



# Specifications

## **General**

Working Length: 65 cm.  
Distal Tip: 7.8 Fr. (2.6 mm)  
Distal Shaft: 9.8 Fr. (3.3 mm)  
Tip Deflection: 160°.  
Working Channel: 3.6 Fr. (1.2 mm)  
Field of View: 80° ± 5.  
Direction of Field: 2-40 mm.

