



Features

- For Arthroscopy and examination/treatment of nasal cavities and nasal pharynx
- Non-Autoclavable; Stainless steel construction to withstand chemical sterilization modalities (EtO, Steris, and Sterrad)
- Wide angle of view and ergonomically designed to directly attach to the camera head for ease of use
- Insertion needle is fully rotatable without the need to reposition the camera head
- For use only with the Smith & Nephew ED3, 360, 460, HD900 and 111 camera heads; Not compatible with the 470 or 560 camera heads



Specifications

Focus Ring	Blue
Focal Length	35 mm
Outer Diameter	2.7 mm
Direction of View	30°
Field Of View	90°
Working Length	67mm
Light Post Direction	Opposite
Cannula Connection	J-Lock
Recommended Light Cable	4.0mm

