



OLYMPUS model CLK-3 endoscopy light source with Air. 150 watt Halogen bulb, small compact unit, easily portable.

Has a dual lamp system allowing instant switch-over, in the event of lamp failure. The built in air pump allows insufflation and lens irrigation and makes it adaptable to all Olympus endoscopes

FEATURES:

- Halogen bulb 150w
- Lamp life approximately 50 hours
- Air feed pump
- Compatible with Olympus OES Endoscopes

SPECIFICATIONS:

- Dimensions: x 7" T x 14 3/4" Long x 8 3/4" Wide by 7" High



Typical Manufacturer's Image