



The Sony UP-51MD is an A5 color video printer using size dye-sublimation technology. The Sony UP-51MD provides approximately 300 dpi (2048 dots) of resolution, 256 gradations for yellow, magenta, and cyan and over 16.7 million colors per dot. The printer features a very fast (20) second print time (Sony UPC-510 media/NTSC mode) and provides RGB, S-Video and Composite Video input and output. Unique to the Sony UP-51MD series printers is Sony's Hi-Eye-Q Quality Imaging System which carefully matches hardware with media products to ensure reliability and consistency.

Features

- 295 DPI with 16.7 Million Colors Per Dot
- 20 Second Print Time
- On-Screen Display of Remaining Ribbon/Print Counts
- Easy Color Setup System
- 8 MB Memory Standard
- Media Choices Include Standard and Self-Laminating



Specifications

Dimensions

Size: 15 5/8 x 5 x 18 3/4

Weight: 31 lbs (approx. 14 kg)

Printing

Paper Size: A5 size - 148mm x 210 mm

Printing Time: NTSC (UP-51MD) - PAL (UP-51MDP)

Memory: 8 Frame Memory

Max paper in tray: 100 sheets

Max. Power consumption: 2.8 A (UC), 1.2 A (CED)