



The Sony UP 870 is Digital black and white printer that is an ideal choice for medical applications including ultrasounds. Includes a UBS interface that allows the user to capture images and then store then on a smaller, portable USB storage device.

Features

- Ideal for today's advance ultrasound systems
- User-friendly controls on the front panel
- Compact, space-saving design
- High resolution of 325 DPI
- High-speed printing of approximately 1.9 seconds
- USB 2.0 Interface
- Utilizes existing 110 Series media by Sony UPP-110S/5, UPP-110HD, UPP-110HG.

