



Sony UP980 (UP-980) video printer Ideal for medical applications including ultrasound, mobile c-arm, and cardiac cath lab, the Sony UP-980 Black and White Video Graphic Printer produces high-quality images in just 10 seconds. A dual scan feature allows the printer to sync up to any 15.75 kHz or 31.5 kHz composite signal. The Sony UP-980 black and white video graphic printer also produces high-quality film-like transparencies with optional blue transparency media.

Features

- True 256 gray scale gradations
- Dual Scan (15.75 or 31.5 khz)
- Multiple pictures on one sheet capability (one, two, four, six)
- Front panel LCD display
- Video loop through
- UL2601-1 listed for medical use
- 10 Second print time



Specifications

Dimensions	W x H x D 316mm x 132.5mm x 305mm
Printer Type	B &W
Printer Gradation	256 gray scale gradation
Power Requirements	100V to 240V, 50/60 Hz
Print Speed	10 seconds / 150mm TSC
Print Driver Software	No
Printing Method	Thermal Graphic Printing

