



Sony UP980 (UP-990) video printer Ideal for medical applications including ultrasound, mobile c-arm, and cardiac cath lab, the Sony UP-990 Black and White Video Graphic Printer produces high-quality images in just 10 seconds. A dual scan feature allows the printer to sync up to any 15.75 kHz or 31.5 kHz composite signal. The Sony UP-990 black and white video graphic printer also produces high-quality film-like transparencies with optional blue transparency media.

Features

- Analog video and digital USB 2.0 interfaces
- Fast print speed - approximately 8 seconds
- Optional blue transparency film (UPT210BL)
- Excellent print quality and contrast
- Front panel LCD display for easy operation
- Multiple images per sheet (two, four, six)
- Fully automatic cutter
- Compact design for easy integration into cart based imaging systems
- Worldwide power supply and video signal compatibility



Specifications

Dimensions	12.4 x 5.2 x 12.0 inches
Printer Type	B &W
Printer Gradation	256 gray scale gradation
Power Requirments	100 to 240V, 50/60 Hz, 1.3 to 2.4A
Print Speed	Approx. 8 second (Std mode)
Print Driver Software	For Windows® XP / 2000 only
Printing Method	Black & White Hybrid Graphic Printer. (Prints film and paper)
Printer Resolution	325 dpi
Safety Compliance	UL60601-1 Listed, CSA C22.2 No. 601.1 GS, GM UL1950 Listed, CSA C22.2 No.950
Warranty	1 Year

