



The Karl Storz 20133020 300 Watt Xenon lamp Cermax xenon lamp by Excelitas is the xenon light source of choice for fiber optic illumination applications such as endoscopy and surgical headlamps. This is a 300 watt xenon arc surgical endoscope medical bulb with a parabolic reflector. The PE300BF light bulb has a 1000 hour life expectancy and a warranty period for up to 500 hours.

Features

- 300 Watt Xenon bulb
- video in/out BNC sockets
- 6-pin Mini-Din-sockets for optional SCB interface



Specifications

Genreal

Watts: 300W

Volts: 14V

Amps: 21.43A

Length: 1.97 inches - 50.0 Millimeters

