



The Stryker 16L is a high flow insufflator that comes with an endoscopic tool that can make exact pressure measurements up to 31 mm Hg along with work with three different gas flow rates. The Stryker 16L insufflator has four different digital displays that can display pressure, current flow rate, and volume of CO<sub>2</sub> consumed.

The Stryker 16L is an endoscopic insufflator. Stryker insufflators are used to inflate a cavity within the body to give more workroom for the surgeons. Most insufflators use CO<sub>2</sub> to inflate the cavity because it is non-flammable, colorless, and dissolves into the body. Insufflators can be used in radiological imaging as well as endoscopic procedures. Soma Tech Intl offers technically and cosmetically refurbished insufflators up to 50% below OEM prices with the same service and warranty as new. Do you have questions about our insufflators? Contact Soma Tech Intl to reach one of our dedicated sales representatives or call 800-GET-SOMA.

