



For enhanced surgical visualization, the Stryker X8000 Light Source complements the Stryker 1188 HD 3-Chip Camera by enhancing its picture with clean, pure light. This maximizes the ability of Stryker's high-definition video equipment to capture the most brilliant and detailed visualization possible for today's surgical procedures.

Specifications:

The X8000's powerful bulb, convenient operations, safety features and unique Align-Lock technology are some of the many outstanding benefits:

- 300-watt xenon elliptical bulb
- Interface with visible LCD screen
- Pure, clean light
- Align-Lock
- Single-handed insertion
- Electronic Scope Sensing Technology

Typical Manufacturer's Image

