



Specifications:

Illumination level setting	Continuous
Mode	Manual, automatic
Power needed, vac	100-120.220-240
Hz	50/60
Dimensions	41x 36 x 16 (16.1x 14.2 x 6.3)
Type	Multipurpose
Connections	
Number/type	1/not specified
Equivalent light cables	Dyonics, Wolf, other adapters available
Video	
Attachment	Yes
Automatic light adjustment	Yes
High-intensity light	Yes
Lamp used	Xenon
Watts	300
# of lamps	1
Color temp, Å° K	6000
Lamp duration, hr	500

