



**Specifications of the Pentax EPM Video Processor**

- **Video Output:** RGBS, Y/C, Composite (NTSC)
- **Light Source:**
  - **Energy:** 300W
  - **Lamp Type:** Xenon
- **Power Supply:**
  - **Voltage:** NTSC 110V
  - **Frequency:** 50-60 Hz
  - **Input Current:** 4.8A
- **Dimensions:**
  - **Width:** 45cm
  - **Height:** 15cm
  - **Length:** 52cm
- **Weight:**
  - **Weight (lbs):** 55.0lbs
  - **Weight (kg):** 25.0Kg

**The average Pentax EPM Endoscopy System has:**

- Pentax EPM-3500 Processor
- Xenon Light Source
- Video Gastroscope
- Video Colonoscope
- High Resolution Monitor



Typical Manufacturer's Image