



The Sony UP21MD (UP-21MD) printer is the ideal choice for the endoscopy and ultrasound markets. This Sony Color Video Printer is very compact and will help to free up space on medical carts. The Sony UP-21MD width at 8 3/8" allows it to be easily integrated, or rack mounted side-by-side with other equipment on medical carts. Although small in size, the Sony UP-21MD printer offers a wide variety of features. Using Sony's Dye Sublimation printing technology, the Sony UP-21MD produces near A6 size prints with approximately 400dpi in approximately 20 seconds, using the Sony UPC-21L media.

Features

- Superb print quality and color reproduction at approximately 400 dpi resolution
- New Compact Size
- Four Frame Memory
- Versatile Analog Video Inputs and Outputs
- Optional remote control capability
- Front Operation



Specifications

Dimensions	Size:8 3/8(W) x 5(H) x 15 5/8(D)
Printer Gadation	100-120 Volts 1.8A, 220-240 Volts 1.0A
Power Requirments	100-120 Volts 1.8A, 220-240 Volts 1.0A
Printer Speed	UPC-21S approx. 16 seconds, UPC-21L approx. 20 seconds
Printer Resolution	325 dpi

