



The Level 1® H-1025 Fast Flow Fluid Warmer provides a rapid flow of warmed fluids, such as crystalloid or blood product, including red blood cells, as volume replacement. Even a small drop in a patient's body temperature can be detrimental – delaying recovery or allowing infection to set in. With Level 1's proven technology, the blood or fluid is inherently free of any "hot spots," as the H-1025 employs a safe, sealed heating system.

## Features

- **Automated Pressure Chambers:** Rigid pressure chambers accommodate standard blood and crystalloid bags, providing a constant 300mmHg pressure for rapid infusion while allowing for fast and easy bag changes
- **Heat Exchanger:** Aluminum heat exchanger transfers heat 1000 times faster than plastic
- **Gas Vent:** The Gas Vent provides the added feature of automated air elimination, which eliminates microbubbles without the need for frequent repriming.
- **Built-In Safety Features:** Monitors reservoir temperature to prevent overheating, with an audible and visual alarm if reservoir temperature reaches 43.9° C.



# Specifications

## General

230 V input voltage  
1400 W power heater  
7 kg weight  
1.4 litre water reservoir volume  
Designed for gravity flow and  
pressure trauma situations

