



Features

- Digital storage of still images, video sequences, and audio files
- Digital alternative to video printer and video recorder
- Compact housing
- Easy and intuitive handling via touch screen, camera head buttons, or footswitch
- Archiving on DVD, CD-ROM, USB stick, or network path, multisession and multipatient
- SDI, S-Video, and composite video inputs
- ALI video signals are transferred directly to the video monitor
- Print-out of still images via HP ink jet printer
- Compatible with Karl Storz Communication Bus (SCB) and OR1 connect series



Specifications

Dimensions

Height: 4" (10.25 cm)
Width: 12" (30.48 cm)
Depth: 16.5" (41.9 cm)
Max Depth (Screen Extended): 25.25" (64.14 cm)

Screen

Length: 8.25" (21 cm)
Width: 6.25" (15.9 cm)