

Maquet Cardiosave



Maquet Datascope CS300



VS



MAQUET
GETINGE GROUP



MAQUET
GETINGE GROUP



Overview

The **Maquet Cardiosave** represents a giant leap forward in functionality and versatility while continuing to deliver the performance and intuitiveness of the datascope pumps you know and trust. With its large touchscreen display, smaller and lighter design, seamless transition from hospital to transport mode, this pump redefines counterpulsation therapy.

Overview

The **Maquet Datascope CS300** has a true one-button start-up, automatically calibrates the fiber-optic pressure sensor, and features improved arrhythmia tracking and response compared to the previous model. The combination of fast balloon pump and superior catheter result in fast inflate/deflate time. The CS300 has the ability to adjust all timing parameters.

Specifications

Maquet Cardiosave

Height: 44 In (111.8 cm)

Width: 22 In (59.9 cm)

Depth: 27 In (68.6 cm)

Weight: 114 lbs (51.8 kg)

Battery: Lithium Ion (qty 2)

Display: LCD 12.1 In Diagonal

Touchscreen: Yes

Maquet Datascope CS300

Height: 43.1 In (109 cm)

Width: 16.8 In (42.7 cm)

Depth: 22.3 In (56.6 cm)

Weight: 185.3 lbs (84.1 kg)

Battery: Sealed Lead-Acid

Display: LCD 10.4 In Diagonal

Touchscreen: No