



The **BFW Maxillume 250-2** has a dual light source with connections for BFW, Wolf, ACMI, Olympus, and Storz with adjustable brightness slider switches. It has 3 Quartz-halogen 250 watt lamps with a color temperature of between 3,200 and 4,000k.

Features

- 250 watts
- 3 Quartz-halogen lamps
- BFW, ACMI, Wolf, Storz, and Olympus light cable connection options.



Specifications

Dimensions

Height: 6.2" (15.75 cm)
Width: 11.7" (29.72 cm)
Depth: 7" (17.78 cm)
Weight: 16.8 lbs (7.6 kg)

Lamps

Type: Quartz-halogen
Watts: 250
of Lamps: 3
Backup Lamps: 1
Color Temperature: 3,200 - 4,000k
Lamp Duration: 50 hours

Power

Power Needed: 100/120 VAC
Hz (Operating): 50/60 Hz

Connections

Type: 2 rotating turrets
Connections: BFW ACMI, Wolf, Storz, Olympus
Configuration: Multipurpose
Mode: Manual

