



The **Olympus PortaView LF-TP** is a lightweight and highly portable intubation fiberscope. It can be used in emergency situations, at the bedside, in an operating room, or in an ICU. It has a large 2.6 mm suction channel that increases the scope's capability to suction to even highly viscous sputa. With a diameter of 5.2 mm, the LF-TP is compatible with single-lumen tracheal tubes with an internal diameter of 6 mm or greater.

Features

- 5.2 mm outer diameter.
- 2.6 mm suction channel.
- Integrated miniature light source.
- Portable.
- High-resolution images with LED technology.
- Compatible with conventional light sources.

