



The Nikon SMZ-1 Stereo colposcope microscope offers a zoom magnification range of 7x-30x with the included 10x/22mm eyepieces. Designed to work as a standard colposcope but it also has a built-in camera for video colposcopy. The camera and video wiring is present but there is no processor or monitor so I am unable to test the video feature.

It is ready to use and comes with the Nikon SMZ-1 Stereo Zoom colposcope head, two 20x/12 eyepieces, 150W Fiberoptic light, green filter, instrument tray and wheeled stand.



Specifications

Eyepiece Nikon Wide Field 10x

Objectives 0.7x - 3x

Total magnification 7x - 30x

Magnification Type Variable zoom

Stand five feet on casters

Illumination Existing light