



The **Gynex Colposcope Single-Step LED**, also referred to as the CO1100-LED, is a device designed for the viewing the tissues of the vagina and cervix by a telescopic system. A colposcope is used to diagnose and examine abnormalities of the vagina and cervix.

Features

- Hand polished optics.
- Adjustable binocular optics.
- Focal length: 320mm
- Field of View: > 22mm diameter
- Depth of Field: > 3mm
- Smooth adjustable controls.
- Working height: 33" - 47".
- Clear, bright, long-lasting LED light source.
- Light intensity 6.5V / 10,000 LUX.

