



The **Neptune 3** waste management system is able to evacuate hazardous fluid and smoke from your surgical site. The Neptune 3 features a vacuum pump that operates at a noise level of 4.3 sones. This means that this is 48% quieter than the Neptune 2. This unit has multiple suction ranges, 7-inch top display, 8.4-inch main display, and can hold an 8-hour ULPA filter for smoke evacuation.

Features

- Suction range between 50-520 mm-Hg.
- 7" top display and 8.4" Main LCD touchscreen display.
- Smoke evacuation capabilities that can an 80-hour ULPA filter with 10 speed settings. The smoke evacuation runs on a separate pump.
- Electronic overflow protection with mechanical backup.

