

Allied Schuco S330A Surgical Aspirators



The **Allied Schuco S330A** surgical aspirator is a reliable tool for newborns and adults. This unit combines lightweight portability with exceptional durability to cater to the diverse needs of both occasional and individual users. It boasts a vibration-free, easily readable regulator gauge, swift-connect tubing for seamless connections, and a hydrophobic bacteria filter for added hygiene. Powered by AC, it ensures a consistently dependable performance.

Features

- High vacuum range for various suction needs – up to 560 mm Hg (22" Hg)
- Vibration-free gauge for easy reading and vacuum control.
- Durable metal canister brackets accommodate a variety of collection containers.
- Collection canisters with improved scratch and flake resistance.
- Flow range over 40 lpm.



Allied Schuco S330A Surgical Aspirators

Specifications

Dimensions

Height: 10 in (25.4 cm)

Width: 8 in (20.3cm)

Length: 8 in (38.1 cm)

Weight: 11.5 lb (5.2 kg)

Operating

Flow Rate: Greater than 40 lpm at regulator inlet

Temperature: 32°F - 104°F (0°C-40°C)

Vacuum Range: Up to 560 mm Hg (22" Hg)

Maximum Vacuum: 560± 3% mm HG

Sound

Sound Level: Less than 58 dBA

UL Approval

UL approved; meets IEC 60601-1

