



The **AMC60B** anesthesia cart is a great solution to support the use of anesthesia in an operating room and help streamline and improve patient care. The AMC60B offers optimal organization, functionality, and spacious drawers. With its durable aluminum construction and 6 drawers, this anesthesia cart allows for an abundance of storage for medicine, supplies, and equipment required for anesthesia administration.

## Features

- Extension table
- 4 hook IV Pole
- Complete aluminum construction
- Drawer dividers
- Seal lock or padlock option
- Castors with brakes
- Comfortable push handle
- Easy glide drawers



## Specifications

**Dimensions**                      **Height:** 44 in (111.8 cm)  
   **Width:** 26.4 in (67.1 cm)  
   **Depth:** 19.2 in (48.7 cm)

**Drawers**                              **6**

**Drawer Heights**                      **3 Drawers:** 4.3 in (10.9 cm)  
   **2 Drawers:** 5.3 in (13.5 cm)  
   **1 Drawer:** 7.8 in (19.8 cm)

**Features**                              Extension table  
   IV Pole  
   Drawer dividers  
   Seal lock or padlock option  
   Castors with brakes  
   Comfortable push handle Easy  
   glide drawers

