

Olympus Berkeley VC-10 Surgical Aspirator



The **Berkeley VC-10** vacuum is a user-friendly, dependable, and low-maintenance system. Crafted for diverse medical environments, this apparatus ensures seamless operation with its subtly integrated canister slots. Engineered with discretion in mind, it aligns perfectly with Berkeley SafeGuard consumables. Its well-thought-out ergonomic structure guarantees effortless configuration, utilization, and portability.

Features

- VC-10 comes as a complete, ready-to-operate system and is designed to accommodate the safety and space requirements of hospital operating room and clinic environments.
- Comes equipped with a dual collection bottle system and a safety trap for automatic vacuum shut off
- The system includes the option for the physician to operate suction with a foot pedal.
- 50% quieter than comparable systems.



Olympus Berkeley VC-10 Surgical Aspirator

Specifications

Dimensions	Width: 18.3 in. (46.4 cm)
	Depth: 25 in. (63.5 cm)
	Height: 38 in. (96.5 cm)
	Weight: 110 lb (49.895 kg)

Operating Voltage	115 Vac, 60 Hz
--------------------------	----------------

Leakage Current	Less than 150 μ a
------------------------	-----------------------

