



The **Baxa MicroFuse Rapid Rate infuser** is intended for intravenous infusion of adenosine. It can deliver a full syringe dose in 6 minutes, and has two infusion rate settings. The MicroFuse infusers are durable and easy to use with simple push-button controls. This infuser operates on two standard alkaline C-cell batteries.

Features

- Compatible with BD, Monoject, Terumo disposable syringes.
- Two infusion rate selections. 6 minute infusions and 8.5 minute infusions.
- Microprocessor controlled constant delivery rate.



Specifications

Dimensions

Height: 4 in (10.2 cm)
Width: 9.6 in (24.4 cm)
Depth: 2.7 in (6.9 cm)
Weight: 1.56 lb (0.7 kg)

Syringe Compatability

Syringe Sizes: 30 to 60 cc
Compatible Syringes: BD, Monoject, and Terumo disposable syringes

Infusion Information

Accuracy of Flow Rate: $\pm 3\%$ over 1.18 in (29.97 mm)
Infusion Rates
Rate 1: 0.533 inches per minute
Rate 2: 0.381 inches per minute
Delivery Times (Full Syringe)
Rate 1: 6 minutes
Rate 2: 8.4 minutes
Occlusion Pressure: < 45 psi (2327 mm Hg) for approved syringes
Back Pressure Effect: None within ± 100 mm Hg pressure

Power

Power Requirement: 3.0 V dc (battery only)
Batteries: Two (2) C-cell, alkaline only

Alarms

Alarms: Occlusion, end of infusion, low battery, return for maintenance and internal malfunction

