

Biodex 056 Series Ultrasound Table



The **Biodex 056 Series Ultrasound Table** is equipped with 5" (13 cm) locking casters for easy mobility. Patient restraints include two body straps, and a powder-coated finish. This manual positioning table offers flexibility with adjustable height, Trendelenburg, and Fowler Back settings, providing optimal comfort for both patients and medical professionals. With a patient load capacity of 350 lb (159 kg), this ultrasound table combines functionality, durability, and user-friendly design for an easier clinical experience. The drop-down foot section accommodates chair positioning or stirrup use, enhancing the table's adaptability. Additionally, an optional IV Pole attachment is available for added convenience.

Features

- Manual positioning for – Height, Trendelenburg & Fowler Back
- Hourglass tabletop narrows to 18" to facilitate probe positioning
- Drop-down foot section for chair positioning or stirrup use.
- IV Pole attachment is "optional"



Biodex 056 Series Ultrasound Table

Specifications

Dimensions

Height Adjustable: 27.5" to 40" (70 to 102 cm)
Width: 24 in (61 cm)
Length: 75 in (191 cm)
Patient Weight Capacity: 350 lb (159 kg)

Mattress

2" thick (5 cm) foam rubber pad. Naugahyde® covering.
Upholstery Color: Standard color is Graphite (Dark grey, #429)

Motions

Height: 29.5" to 37" (75 to 94 cm), hydraulic foot pump,
Trendelenburg: 0° to 15°, hand crank
Fowler Back: 0° to 80° in five preset positions, manual

Leg Support

Drop Down Leg Section: 14.5" | (37 cm)
Ankle Straps

