

Biodex Echo Pro Echocardiography Table



The **Biodex Echo Pro Echocardiography Table** reduces the total amount of time required to achieve a quality image. Patient comfort is paramount, with cushioned arm supports and a thick mattress ensuring a comfortable experience throughout lengthy procedures. The cardiac scanning cutout provides open access to the patient's left thorax area for an unobstructed apical approach. Encounter a smooth interaction between sonographer and patient, fostering correct scanning posture to mitigate musculoskeletal injuries.

Features

- Allows seated or standing work positions
- Optional vascular scanning arm board mounts to either side of table
- Adjustable sonographer cutout for maximum patient access and comfortable scanning positioning
- Motorized height, Fowler positioning
- Swing-out cardiac scanning cutout doubles as armrest



Biodex Echo Pro Echocardiography Table

Specifications

Dimensions

Height Adjustable: 23 in to 39 in (58.4 to 99 cm)
Width: 30 in (76.2 cm)
Length: 73 in (185.4 cm)
Weight: 335 lb (151.9 kg)
Patient Weight Capacity: 500lb (226.7 kg)

Mattress

Torso section 2" (5 cm) thick, center and leg sections 3" (8 cm) thick; foam rubber pad covered with seamless, durable Naugahyde®

Upholstery Color: Graphite

Motions

Height Adjustable: 23" to 39" (58.4 to 99 cm)
Trendelenburg: 0° to ±15°
Fowler Back: 0° to 80° infinitely adjustable

Controls

Hand Control: Activates height, Trendelenburg motions, Fowler positioning and auto level motions
Foot Control (optional): Activates height, Trendelenburg motions, and Fowler positioning

Wheels

5" (12.7 cm) Individual locking swivel casters; central floor locking system

