

Biodex Ultra Pro Ultrasound Table



The **Biodex Ultra Pro Ultrasound Table** is a specialized table for all general ultrasound procedures. With an ergonomic design catering to both patients and sonographers, this table enhances the overall efficiency of achieving high-quality images. The open frame design facilitates closer proximity for sonographers, particularly beneficial for tasks like carotid imaging or long leg arterial scanning. Equipped with a motorized Fowler back, this table offers infinite adjustability up to 80° through a hand control or an optional foot controller.

Features

- The Central Floor Locking System secures table position and is accessible from either side of the table
- Foot Support
- Adjustable sonographer cutout for maximum patient access and comfortable scanning positioning
- Motorized height, Fowler positioning, Trendelenburg motions and auto-level adjustment



Biodex Ultra Pro Ultrasound Table

Specifications

Dimensions

Height Adjustable: 23 in to 39 in (58.4 to 99 cm)
Width: 30 in (76.2 cm)
Length: 70 in (177.8 cm)
Weight: 335 lb (151.9 kg)
Patient Weight Capacity: 500lb (226.7 kg)

Mattress

Torso section 2" (5 cm) thick, center and leg sections 3" (8 cm) thick; foam rubber pad covered with seamless, durable Naugahyde®

Upholstery Color: Graphite

Motions

Height Adjustable: 23 in to 39 in (58.4 to 99 cm)
Trendelenburg: 0° to ±15°
Fowler Back: 0° to 80° infinitely adjustable

Controls

Hand Control: Activates height, Trendelenburg motions, Fowler positioning and auto level motions
Foot Control (optional): Activates height, Trendelenburg motions, and Fowler positioning

Wheels

5" (12.7 cm) Individual locking swivel casters; central floor locking system

