



The BOVIE Specialist is a versatile and safe unit providing a unique coagulation/hemostasis waveshape with superior clinical performance. The availability of accessories and a wide selection of electrodes will provide many surgeons with an instrument that can be utilized effectively in specialized minor surgery. The Specialist is designed with a graphic around the Coagulating Power Control Knob which depicts the power settings to be used for various types of electrodes. The graphic is based on bi-terminal techniques utilizing a dispersive electrode. It is not an exact setting for any procedure. However, it can be used as an aid for safe approximate settings, thereby reducing the chance of electrode damage or excessive tissue destruction.

Features

- Circuit Sentry monitoring (Dispersive Electrode).
- Coagulation overrides when inadvertent actuation of cut and coagulation is simultaneous.
- Coagulation electrode symbols.
- Output only to the selected active electrode.
- Distinctive separate audio tones.
- Bipolar outlet.



Specifications

Dimensions

Height: 4.75 in (12 cm)
Width: 10 in (25.4 cm)
Depth: 14.25 in (36 cm)
Weight: 20 lbs. (9 kg)

Monopolar Power Consumption

Cut: 75 W
Blend: 50 W
Coagulation: 36 W

Bipolar Power Consumption

Cut: 32 W
Coagulation: 36 W

Operating Frequency:

1.75 MHz

Line Voltage:

90 - 130 V

Line Frequency:

50 - 60 Hz

Current:

3.2 A

