

Covidien Kangaroo Connect Feeding Pump



The **Covidien Kangaroo Connect Feeding Pump** is a portable solution for delivering enteral nutrition with precision and flexibility. Suitable for various medical environments including hospitals, acute care facilities, and home care settings, it ensures reliable performance in both stationary and ambulatory conditions. It offers continuous, dose, or bolus feeding options when used with Kangaroo Connect feeding sets. Compatible with essential accessories and optional add-ons, it ensures reliable performance for patients from infants to adults, excluding neonates.

Features

- Quiet pump operation
- Large, color LCD display
- Compact, lightweight design
- Rubberized casing for better grip
- View previous 72 hours of feeding history
- Pole clamp can detach from either pump or WCH without removing from IV pole



Covidien Kangaroo Connect Feeding Pump

Specifications

Dimensions

Pump

Height: 3.9 in (9.9 cm)

Width: 6.1 in (15.4 cm)

Depth: 1.6 in (4 cm)

WCH

Height: 5 in (2.8 cm)

Width: 7.6 in (19.2 cm)

Depth: 3.5 in (9 cm)

Weight: 73 lbs (0.33 kg), 1.4 lbs (0.635 kg) with pole clamp

Material

Soft-Touch Coating: Latex-Free Thermoplastic Urethane Pump Housing: Flame resistant Polyester/Polycarbonate blend WCH Housing: Flame resistant ABS/Polycarbonate blend

Occlusion Pressure

Maximum Occlusion Pressure: 20 psi (138 kPa)

Accuracy

5% or 0.5 mL/hr, whichever is larger, for all delivery rates no matter the type of Kangaroo Connect feeding set. The top of the fluid column should be at a starting height of 25.4 cm (10") \pm 0.76 cm (0.3") above the top of the pump.

Feeding Formula Dose

1 - 3000 mL in 1 mL increments

Pumping Mechanism

Rotary Peristaltic

