

# Medtronic Covidien WarmTouch WT 6000 Patient Warmer



The **Medtronic Covidien WarmTouch WT 6000** patient warmer is a convective warming unit equipped with an intuitive screen for monitoring temperature across five settings. Engineered to maintain normothermia effortlessly, this system ensures patient comfort and optimal airflow with its soft yet durable WarmTouch blanket, ideal for various surgical procedures. Experience easier temperature management and enhanced patient care with the WarmTouch WT 6000.

## Features

- LCD screen
- Onscreen filter expiration warning
- Lightweight, easy to transport and attach the unit to a cart, bed rail, or an IV pole.



# Medtronic Covidien WarmTouch WT 6000 Patient Warmer

## Specifications

### Dimensions

#### Warming Unit -

**Height:** 24 in (60 cm)

**Width:** 17 in (43 cm)

**Depth:** 12 in (30 cm)

**Weight with power cord:** 11.5 lb (5.2 kg)

#### Transport Cart (Optional) -

**Height:** 26.4 in (67.1 cm)

**Width:** 12.7 in (32.3 cm)

**Depth:** 15.2 in (38.6 cm)

**Weight:** 6.8 lb (3.1 kg)

### Maximum contact surface temperature

44.1°C (111.4°F)

### Thermal protection threshold

Thermostat (internal) 49°C to 55°C (120°F to 131°F)

### Power requirements

#### 100 to 240 V AC

Max. current at 100 V = 8 A

Max. current at 240 V = 5 A

### Temperature

18°C to 28°C (64.4°F to 82.4°F)

### Input frequency

50/60 Hz ±1 Hz

### Average alarm level

56 dB @ 1 meter

### Automatic temperature step down

After 45 minutes of continuous use, blower will step down from the boost to high setting.

