

# Midmark Ritter 646/647 Podiatry Chair Procedure and Exam Tables



The **Midmark Ritter 646/647 Podiatry Chair** is designed to meet the demands of lengthy procedures. Its ergonomic foot section provides physicians with essential support at the point of care. With easy release handles and integrated touch controls, this chair streamlines your workflow, while custom-designed motors swiftly position patients for maximum productivity and ergonomic comfort.

## Features

- Powered base, back and tilt for effortless positioning. Easy-to-reach release handles on both sides of the foot section allow for smooth, one-handed operation.
- Large stainless steel debris tray.
- Adjustable chair arms rotate completely out of the way for improved comfort and accessibility.



# Midmark Ritter 646/647 Podiatry Chair Procedure and Exam Tables

## Specifications

### Dimensions

#### Model 646

- Base Height: Fixed height of 21" (53.3 cm)
- Foot Section Height: With foot retracted 15.5" to 54" (39.4 cm to 137.2 cm)
- Seat Width with Chair Arms Installed: 37.5" (95.3 cm)
- Patient Weight Capacity: 450 lb (204 kg)

#### Model 647

- Base Height: Minimum 19" (48.3 cm) Maximum 40" (101.6 cm)
- Foot Section Height: With foot retracted 15.5" to 54" (39.4 cm to 137.2 cm)
- Seat Width with Chair Arms Installed: 37.5" (95.3 cm)
- Patient Weight Capacity: 450 lb (204 kg)

### Rotation

**Back Support Angle Range:** 0° to 80

**Seat Tilt Angle Range:** 0° to 30°

Allows 350° of chair rotation and adds 3.5" to the seat height specifications

### Electrical

115 VAC, 6.0 A, 60 Hz

