



The **Blickman 7925** is a tabletop Blanket warmer with 2 variants: A glass door and a solid steel door. This warming cabinet can sit on an existing countertop, or on a mobile stand. It comes with an adjustable perforated stainless steel shelf. It is primarily used for warming solutions.

FEATURES

- Comes in two variants: Glass door, and Steel door.
- The unit can be programmed to heat only between 90°F - 110°F
- Can be used as a tabletop, or under-counter warming cabinet
- Easy-to-use keypad with controls for setting internal temperature



Specifications

Dimensions

Height: 30.0 in (76.2 cm)
Width: 24.5 in (62.2 cm)
Depth: 26.6 in (67.5 cm)
Weight: 215 lbs. (97.5 kg)
Internal Cubic Foot Capacity: 5.05
Shelves: 1

Electrical

Total Wattage: 750
Amps: 7.5

Materials

[304 Stainless Steel](#)
[1/4" clear double panel, insulated, tempered sheet glass](#)
[Perforated Stainelss Steel Shelf](#)

