

# Conmed Beamer CE200 Electrosurgical Unit



The **Conmed Beamer CE200 Electrosurgical Unit** merges advanced technology with user-friendly design. Engineered to accommodate a wide array of endoscopic therapies within a singular platform, this system empowers practitioners with significant versatility.

## Features

- Hands-Free Double Foot Pedal
- Convenient Storage
- Customizable User Settings
- Easy to use Interface



## Specifications

### Dimensions

**Height:** 7.16 in (18.2 cm)

**Width:** 15.35 in (39 cm)

**Depth:** 17.12595 in (43.5 cm)

**Weight:** 18.29 lbs (8.3 kg)

### LF Leakage Currents

Complying to IEC 60601-1/IEC 60601-2-2, Part 1 (test in connection with Beamer System)

### Duty Type

INT 10 s / 30 s ( = duty factor of 25% )

### Approval / Mark of Conformity

CE 0297, in compliance with 93/42/EEC

### Output Data: Cutting Power (Pulse)

Max. 350 watts

### Output Data: Coagulation Power

Max. 120 watts

### Nominal Frequencies (Mains)

50 / 60 Hz

### Power Req / Mains Current / Fuse

100 – 240 V +/- 10% / max. 6.3 A / 6.3 A

