



The **Medline BioCon 700** is a Bladder Scanner, which utilizes ultrasound technology to accurately assess bladder volume. The device is portable and user-friendly, allowing for quick and efficient bladder volume measurements using BladderPoint technology without the need for catheterization. It helps healthcare professionals assess urinary retention, monitor urinary output, and manage bladder-related conditions, ultimately improving patient care and comfort.

## Features

- Quickly and confidently locate the bladder for accurate measurement with BladderPoint technology that puts a bull's eye in the middle of the bladder
- Dependable performance with built-in self-diagnostics to alert you of any hardware or maintenance issues
- Tough polycarbonate casing promotes consistent, uninterrupted service and is supported by a five-year warranty including drop protection
- Easy to use on bariatric patients with rear grip on the probe and no need to fan probe



# Specifications

**Dimensions**

**Height:** 4.1 in (10.3 cm)  
**Width:** 7.95 in (20.2 cm)  
**Depth:** 2.6 in (36.6 cm)  
**Weight:** 2.0 lbs. (0.92 kg)

**Power**

**Power Supply:** 3V, 9V DC adapter  
**Power Input:** AC100-240V 50/60Hz  
**Battery:** Integrated 7.4 Li-ion rechargeable battery

**Other**

**Display:** 7" TFT-LCD (WVGA); 16.7M colors  
**Range/accuracy:** 0-999 ml;  $\pm 15\%$ ,  $\pm 15\text{ml}$   
**Storage capacity:** 700 scans

