

Portable Aspirator

# ALLIED GOMCO 4005

## ABOUT

The Gomco 4005 Aspirator Pump is a highly efficient, general-purpose suction device suitable for various medical applications, including emergency rooms, recovery rooms, outpatient clinics, and dental practices.

## FEATURES

**Suction System:** Equipped with a gauge and a suction regulating valve, this unit can generate up to 25 inches Hg of vacuum.

**Safety and Compliance:** The unit is listed as UL/ULC, guaranteeing adherence to safety standards. It also comes with a hospital-grade plug.

**Large Capacity Canister:** The Gomco 4005 features a graduated 11500 ml disposable canister, allowing for extended use without frequent changes.

**Contamination Prevention:** To enhance safety, the unit includes a filter that helps prevent fluid and aerosol contamination.



# SPECIFICATIONS

## DIMENSIONS

---

**Height:** 9.5 in (24.13 cm)

**Width:** 12.5 in (31.75 cm)

**Depth:** 10.0 in (25.4 cm)

**Weight:** 14.5 lb (7 kg)

## POWER

---

**Voltage:** 115 V / 60 Hz

**Electrical Cord:** 6 in Hospital grade

## GENERAL

---

**Vacuum Range:** -80 to -635 mm Hg

**Flow Rate:** 40 lpm at open flow

**Sound Level:** Less than 54 dBA

**Pump Type:** Carbon rotary vane

**Disposable Canister:** 1 - 1500 ml