

BAXTER®/BARD® INFUSO.R™ PUMP

- Designed to aid in many intravenous agents given during anesthetic procedures.
- Provides for the convenient delivery of narcotics, muscle relaxants, vasoactives and other drugs routinely given during anesthesia
- Easy to use and set up
- Drug delivery choices are easily selected and changed using front panel rotary switches



Typical Manufacturer's Picture

Specifications

Accuracy	Infusion: linear rate $\pm 3\%$. Bolus: linear displacement $\pm 3\%$.
Battery Life	Typical use: 150 hours. At low battery and typical use: 4 hours
Batteries	Four "C" – size Alkaline Cells
Operating Voltage	5 – 7 Volts
Operating Current	20 mA
Low Battery Voltage	5.0 Volts approximately
Flow Rates	0 to 600 mL/hr
Flow Profile	Bolus: continuous. Infusion: pulsed intermittently
Occlusion Force	8 ± 1 lbs.
Maximum Occlusion Pressure	11 psi w/60 cc syringe; 20 psi w/20 cc syringe
Occlusion Detection Time	120 seconds at 36 mL/hr (inverse related to flow rate)
Volume Stored on Occlusion	1.1 mL approximately
Back Pressure	
Effect on Accuracy	None to Occlusion Pressure
Size	9.2 x 4.5 x 2.0 inches
Weight	2 lbs. With batteries
Syringes	Monoject® or B_D® 60 cc; Monoject® or B-D® 20 cc

Soma Technology, Inc. acknowledges all registered trademarks of manufacturers' listed.
The technical data given in this publication are for general information and are subject to change without notice.



166 Highland Park Dr. Bloomfield, CT 06002 U.S.A.
Tel: (860) 218-2575 Fax: (860) 218-2565 Email: soma@somatechnology.com