



The **GE 6Tc-RS** is a transesophageal echocardiography (TEE) probe that is compatible with the GE Vivid i system and designed to be used for adult cardiac applications. If the high cost of a new 6Tc-RS TEE probe is preventing you from purchasing one, then you need to check the price of a refurbished 6Tc-RS probe at SomaTechnology.com.

Features

- **Bandwidth:** 2.9 – 8 MHz
- **Applications:** Cardiac
- **Imaging Modes:** 2D, Harmonics, M-Mode, Color Doppler, PW Doppler, CW Doppler, Contrast Imaging, TVI, Tissue Tracking, TSI
- **Footprint:** L45 x W14 (x H12.5) mm (Logiq e), 12 x 14 mm, length 45 mm (Vivid i, q)
- **Field of View:** 90°
- **Compatibility:** Logiq e, Vivid E9, Vivid i, Vivid i n, Vivid q

