



The **McGRATH Series 5** video laryngoscope is a device that aids in the intubation of the trachea. A trachea blade restrains the patient's tongue, which allows the clinician to have a clear view of the larynx and entrance to the trachea. It features a light source and small camera within the CameraStick assembly that allows full view of the larynx during an intubation procedure. One AA 1.5V battery mounted in the handle of the device powers the LCD screen, camera, and light source.

Features

- Adjustable blade.
- Disposable intubation blades.
- 1.7" Color LCD Display.
- Weight: 0.6 lb



Specifications

Dimensions (Laryngoscope Assembly)

Height: 7.9" (20 cm)
Width: 2" (5 cm)
Depth: 5.9" (15 cm)
Weight: 0.6 lbs (0.27 kg)

Disposable Intubation Blade

Height: 4.3" (11 cm)
Width: 0.5" (1.2 cm)
Depth: 0.6" (1.5 cm)
Material: Medical grade optical polymer

General

Power: AA 1.5V battery giving typically 1 hour of non-continuous use.
Light Source: High intensity LED
Display: Color LCD (1.7" diagonal)
Camera: VGA
Protection: IP65
Materials: Durable medical grade rubber, thermoplastics and stainless steel. This device and packaging are latex free.

