



The **Alaris CareFusion 8300 EtCo₂ Module** provides continuous respiratory monitoring. The Infusion Pump monitors CO₂, EtCO₂ (exhalation) FiCO₂ (inhalation), and respiratory rate. The Carefusion EtCO₂ Module is a side-stream capnography. The Disposables deliver a sample of the inhaled and exhaled gasses from the ventilator disposable or directly from the patient into the monitor for CO₂ measurement moisture and patient secretions are extracted from the sample by the micro stream inline filter while maintaining the shape of the CO₂ waveform. The Alaris CareFusion 8300 EtCO₂ Module release a gas sample through a micro-sample cell (15 Microliters).

Features

- Reduce risks of opioid infusions
- Customizable profiles add flexibility across all patient care areas
- Utilizes sampling lines that cuts the need for moving parts, external sensors or routine calibration
- Provides patient EtCO₂ sampling products to meet a broad spectrum of sampling needs available with or without oxygen delivery

