



The **Alaris IMED Gemini PC-1** is a single channel volumetric infusion pump and controller which provides accurate and automatic infusion of intravenous drugs and fluids. The Gemini PC-1 utilizes peristaltic action in both the pump and controller modes. When in controller mode, it senses the pressure created by the height of the fluid container and uses that pressure as the occlusion pressure limit. If in pump mode, the occlusion pressure limit is preset at about 10 ( $\pm 2$ ) psi. The infusion rate is selectable from 0.1 mL/hr to 999 mL/hr.

## Features

- Can be operated as a pump or controller.
- Secondary (piggyback) infusion can be administered using infusion parameters independent of the primary infusion.
- Infusion rate selectable between 0.1 mL/hr and 999 mL/hr.
- Large easy-to-read display.



