

Features

-
- A blue manual suture driver by Allen Medical Systems. It features a long, thin metal shaft with a black handle at the bottom and a blue handle at the top. The top handle has a trigger mechanism and a small label that reads "Allen Medical Systems". The bottom handle has a black trigger and a small label that reads "Allen Medical Systems". The device is shown against a white background.

Specifications

General

Patient Weight Capacity: 350 lbs (159 kg)

Lithotomy Range: +90° to -22°

Abduction Range: +25° to -9°

Procedures

GYN, Urological, Laparoscopic, Robotic, and Colorectal procedures including:

TAVH, LAVH

Colposcopy

Cryosurgery

LEEP Procedure

D&C

Radical Perineal Prostatectomy

General Surgery

AP Resections

Cysto Procedures

Radical Perineal Prostatectomy

Oncology Procedures

Robotic Surgical Procedures

