



The **Allen Ultrafins** stirrups is an upgrade of the Yellofins stirrup. It has a reinforced rod and joint assembly and the flexible boot that helps the whole system be able to handle a patient weighing up to 800 lbs (363 kg). The Allen Ultrafins are still appropriate for use with general patients even though the Ultrafins were designed for large patients.

Features

- Flexible boot that comfortably contours to the calf.
- Unique fin design that covers the fibular head to protect superficial peroneal nerve.
- Squeeze grip handle that provides easy intraoperative adjustment.
- Lithotomy and length Indicators.



Specifications

General

Patient Weight Capacity: 800 lbs (363 kg)

Lithotomy Range: +84° to -33°

Abduction Range: +25° to -9°

Procedures

GYN, Urological, Laproscopic, Robotic, and Colorectal procedures including:

TAVH, LAVH

Colposcopy

Cryosurgery

LEEP Procedure

D&C

Radical Perineal Prostatectomy

General Surgery

AP Resections

Cysto Procedures

Radical Perineal Prostatectomy

Oncology Procedures

Robotic Surgical Procedures

